Introduction to the Gait Trainer

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All tests, if not otherwise specified, were performed at an ambient temperature of $21\text{ }^{\circ}\text{C} \pm 5\text{ }^{\circ}\text{C}$.

All tests were performed with the height adjustments at their maximum and swivelling wheels in the least stable position. The angle was always recorded.

Parking or running brakes shall not be activated unless specified in the test procedure.

During the stability tests, the wheels were in the most unstable position for the test being performed. The gait trainers were prevented from sliding or rolling before tilting occurred. The results of the tests were not influenced by the means used.

All measurements are outside dimensions except for the wheel width

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Stability Tests

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Forward

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Place the gait trainer with its wheels on a plane that can be tilted from the horizontal so that the gait trainer is facing downwards when the plane is tilted. The centre line of the hinges were parallel to the axis of the front wheels. The front swivelling wheels were positioned with the ground contact point to the rear of their castor axis.

The position of the load was $60\text{ mm} \pm 5\text{ mm}$ behind the front end of the gait trainer. The load must be applied vertically and pass through the midpoint of the supports. The load must weigh $250\text{ N} \pm 2\%$.

Tilt the plane and record the maximum angle of the plane at the point of tilting, rounded down to the nearest 0.1° .

The angle of the plane at the point of tilting must exceed 15.0° .

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Backward Stability

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Place the gait trainer with its wheels on a plane that can be tilted from the horizontal so that the gait trainer is facing backwards when the plane is tilted. The centre line of the hinges must be parallel to the axis of the front wheels. The front swivelling wheels were positioned with the ground contact point to the rear of their castor axis.

The position of the load must be $60\text{ mm} \pm 5\text{ mm}$ in front of the back of the gait trainer. The load must be applied vertically and pass through the midpoint of the supports. The load must weigh $250\text{ N} \pm 2\%$.

Tilt the plane and record the maximum angle at the point of tilting, rounded down to the nearest 0.1° .

The angle of the plane at the point of tilting must exceed 7.0° .

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Lateral Stability

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Place the gait trainer with its wheels on a plane that can be tilted from the horizontal so that the gait trainer faces across the plane. The centreline of the hinges must be parallel to the line through the centres of the wheels on the side of the gait trainer which is closest to the hinges. The position of the load must pass through the centre of that forearm support that is closest to the hinge. Apply the load vertically. The load must weigh $250\text{ N} \pm 2\%$.

Tilt the plane and record the maximum angle at the point of tilting, rounded down to the nearest 0.1° .

The angle of the plane at the point of tilting must exceed 4.5° .

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Results of the stability testing can be found [here](#)

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Strength Tests

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The Gait Trainer must be placed with its wheels on a horizontal and stationary surface.

The load must be applied vertically to the forearm supports of the Gait Trainers. The force vector must be vertical and pass through the middle of the forearm supports.

A load was gradually applied. The maximum load of **76 kg (167 lbs)** was applied for **10 minutes**.

Critical measurements were documented pre-loading, during loading and after the load was removed. The Gait Trainers must be examined for any cracks, fractures and potential hazards.

Results of the strength testing can be found [here](#)

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The stability testing on the gait trainers was carried out in accordance with ISO 11199-3 (version 2005): "Walking aids manipulated by both arms -Requirements and test methods (Part 3 - Walking tables)"

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The strength testing on the gait trainers was carried out in accordance with ISO 11199-3 (version 2005): "Walking aids manipulated by both arms -Requirements and test methods (Part 3 - Walking tables)"

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*Starting from the front left leg, then each leg in turn (clockwise)

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